

Blood Pressure Monitor

SPECIFICATION

- Measurement Range -
- Pressure - 0mmHg & 500mmHg
- Pulse - 30 to 180 Beats/Minutes
- Measurement Accuracy Pressure - 0.3mmHg
- Pulse - 0.5%
- Pressurization- Automatic Pressurization
- Memory - 120 memories in four Groups with date and time
- Function - Low Battery Detection Automatic Power-off
- Power Source - 4 AA Batteries
- Unit Weight -
- Approx 405g(14.30o.z)
- (Excluding Battery)
- Unit Dimensions - Approz 110 x 115 x 70mm (4.33" x 6.10" x 2.76")(L x W x H)

